

3M™ ESPE™ Paradigm™ MZ100 Block for CEREC®

New independent studies
prove the benefits of using

3M™ ESPE™ Paradigm™ MZ100 Block for
CEREC®, including its enamel-like wear
and high fracture toughness.



Studies presented at the
2000 and 2001 IADR/AADR
meetings in Chicago and
Washington compared the
bond strengths, fracture
toughness, clinical
performance, and wear
on enamel. The following
graphs illustrate some of
the results. Study results
have been published in the
Journal of Dental Research,
2000, v. 79, Special Issue:
Abstracts of Papers, and
Journal of Dental Research,
2001, v. 80, Special Issue:

Abstracts of
Papers.




YF Mansour, CA Mitchell

Division of Restorative Dentistry, Queen's University, Belfast, U.K.: Effect of Thermocycling on Fracture Toughness of Porcelain and Polyglass Materials. J Dent Res 80 (AADR Abstracts #570) 2001.

This study measured the effect of thermocycling on the fracture toughness (K_{IC}) of dental porcelains and polyglasses used to fabricate bonded aesthetic restorations, using the chevron notch short rod method. Samples were stored isothermally (37°C, 7 days) or thermocycled (10,000 cycles, 5°C–60°C). Thermocycling had an statistically significant effect on toughness. In addition, differences were seen between some of the materials tested.

FRACTURE TOUGHNESS




K-H Kunzelmann, J Manhart, A Mehl, R Hickel

Department of Restorative Dentistry, LMU University, Munich, Germany: The Effect of CEREC CAD/CAM Materials on the Wear of Human Enamel. J Dent Res 80 (AADR Abstracts #575) 2001.

The wear of human enamel against composite and porcelain materials for CEREC was measured on a newly developed wear simulator. The wear test comprised 50,000 cycles of a vertical load of 50 N in distilled water at room temperature. Laser scans of silicone replicas were used to measure the wear of the enamel and materials.

"The total wear (material+enamel) of the composite material is far less than the total wear of the ceramic materials. No significant differences were found among the CEREC and lab-made ceramics."

TOTAL VOLUME LOSS



3M ESPE

3M™ ESPE™ Paradigm™ MZ100 Block for CEREC®


DJ Fasbinder, JB Dennison, D Heys, K Lampe

University of Michigan, Ann Arbor, MI, and University of Zurich, Zurich, Switzerland: Clinical Evaluation of CAD/CIM-Generated Polymer Ceramic Inlays. J Dent Res 80 (AADR Abstracts #1882) 2001.

Six-month results are presented of a longitudinal, randomized clinical trial comparing Paradigm MZ100 block and Vita™ Mark II inlays. Restorations of both materials were fabricated on a CEREC 2 and bonded with 3M™ ESPE™ RelyX™ ARC Cement. In response to cold stimulus, no sensitivity was reported for Paradigm MZ100 block at one week, and one case for Vita™ Mark II. Restorations were evaluated with a modified USPHS rating at baseline and six months.

"No sensitivity was reported in either group after two weeks, and six months. All restorations were rated clinically acceptable at six months."

CLINICAL EVALUATION




JR Sturdevant, EJ Swift Jr., SC Bayne

University of North Carolina, Chapel Hill, NC: Cement Bond Strength to Millable Composite for CAD/CAM Restorations. J Dent Res 80 (IADR Abstracts #2479) 2000.

In this study, the adhesion of RelyX ARC cement to Paradigm MZ100 machinable composite block is measured with respect to variables: roughening with air abrasion vs. sanded, and ethyl acetate rinse vs. no solvent wash. The 180 grit sanding mimics the standard as-milled surface. Buttons of mill block were cemented to mill block substrate with these treatments; all surfaces were treated with 3M™ ESPE™ RelyX™ Ceramic Primer, per product instructions. The results show that the solvent had no effect on adhesion. While the adhesion with the sanded surface was excellent, air abrasion notably improved the adhesion.

BOND STRENGTH



3M ESPE Technical Hotline: 1-800-634-2249

3M ESPE Web Site: www.3MESPE.com

3M, Paradigm and RelyX are trademarks of 3M.

ESPE is a trademark of 3M ESPE AG.

CEREC is a registered trademark of Sirona Dental Systems GmbH.

Vita is a trademark of Vita Zahnfabrik H. Rauter GmbH & Co.

3M ESPE

Dental Products

3M Center
Building 275-2SE-03
St. Paul, MN 55144-1000
USA

Dental Products

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1 800 265-1840 ext. 6229



40% Pre-consumer waste paper
10% Post-consumer waste paper

Printed in U.S.A.
Copyright © 3M (IPC) 2001.
All rights reserved.
70-2009-2352-5