Instructions for Use Mode d'emploi Instrucciones de uso

3M ESPE AG Dental Products D-82229 Seefeld – Germany
3M ESPE Dental Products St. Paul, MN 55144-1000 U.S.A.
3M ESPE Customer Care/MSDS Information in U.S.A. 1-800-634-2249 and Canada 1-888-363-3685.
3M, ESPE, Filtek, Garant, Penta, Protemp, RelyX and Sof-Lex are trademarks of 3M or 3M ESPE AG. © 3M 2010. All rights reserved.

en	ENGLISH
fr	
es	

Product Description
Protemp™ Plus, manufactured by 3M ESPE, is a composite material for the fabrication of temporary restoration at the chairside. Based on multi-functional methacrylate esters, this two-component system has filling material-like properties and is available in the following colors: A1; A2; A3; A3.5*; B3; Bleach. With the Garant™ Dispenser, manufactured for 3M ESPE, the material can be applied directly from the cartridge without causing any bubbles to form. Protemp™ Plus temporary restorations can be modified and individualized using low-viscosity composites (e.g., from the Filtek™ Flowable Composite product line, manufactured by 3M ESPE). * Not available in all countries.
For details on all mentioned products, please refer to the corresponding Instructions for Use. Instructions for Use should not be discarded for the duration of product use.

Application Areas
- Fabrication of temporary crowns, bridges, inlays, onlays, and veneers - Fabrication of long-lasting temporary restorations - Lining material for prefabricated temporary crowns made of composite (e.g., Protemp™ Crown, manufactured by 3M ESPE) and metal (e.g., Iso-Form crowns, manufactured by 3M ESPE)

Precautionary Measures
For Patients
- Avoid use of this product in patients with known acrylate allergies. - This product contains substances that may cause an allergic reaction by skin contact in certain individuals. - If prolonged contact with oral soft tissue occurs, flush with large amounts of water. - If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

For Dental Personnel
- This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. - If skin contact occurs, wash skin with soap and water. - Use of protective gloves and a no-touch technique is recommended. - Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If contact with eyes or prolonged contact with mucous membranes occurs, rinse thoroughly with plenty of water. - If symptoms persist, seek medical attention.

MSDSs for the mentioned 3M ESPE products can be obtained at www.mmm.com or from your local subsidiary.

Preparatory Work
- Take an alginate, silicone, or polyether impression from the intact row of teeth (e.g., using a product manufactured by 3M ESPE). - Instead of an impression, a laboratory-made vacuum-formed template or, for single crowns, a preformed crown can be used for shaping. - To enhance the strength of the temporary restoration, remove interproximal grooves in the impression. If required due to tight spacing conditions, enlarge the relevant sites by cutting. - For easy replacement in the mouth, relieve undercuts. - With missing teeth or gaps in the molar area, cut a groove in the impression to produce a stable, bar-shaped connection in the temporary restoration. - If required, prior to taking the impression close any gaps between front teeth with denture teeth acting as spacers; interlock multiple denture teeth with wax.

Times

Option 1 for optimum surface smoothness and gloss

00:00–00:40 min:sec	Place Protemp™ Plus in the impression and reposition in the mouth
00:40–01:40 min:sec	Setting in the mouth
01:40–02:50 min:sec	Best point in time for removal of the impression
Up to 05:00 min:sec	Complete setting
From 05:00 min:sec	Remove the temporary restoration from the matrix. Remove oxygen inhibition layer. Final finishing

Option 2 for faster results

00:00–00:40 min:sec	Place Protemp™ Plus in the impression and reposition in the mouth
00:40–01:40 min:sec	Setting in the mouth
01:40–02:50 min:sec	Remove from mouth and matrix, trim margins
From 05:00 min:sec	Remove oxygen inhibition layer. Final finishing

This timeline will deliver faster results. For optimum surface smoothness and gloss, follow the timeline option 1.

The processing times are applicable to products which have been stored and processed at a temperature of 23 °C/74 °F and at 50% relative humidity. Higher temperatures shorten these times while lower temperatures extend them.

Dosing and Mixing
The pastes are dosed and statically mixed in the Garant dispenser. Use the enclosed Garant™ Mixing tips blue only! Garant Mixing tips are manu-factured for 3M ESPE.
During the Initial Dosing
Remove and discard the cartridge lid. Under no circumstances should the cartridge lid be used to reseal the cartridge as this could result in the inclusion of air bubbles. - Check whether both cartridge openings are free from obstruction and remove obstructing material, if any, with a suitable instrument. - Before a mixing tip is applied for the first time, squeeze paste out onto a mixing block until the base and catalyst pastes are the same. Discard the paste which has been squeezed out.

Application
It takes the material a certain amount of time to emerge from the newly attached mixing tip. The flow of material is interrupted as soon as the pressure on the

handle is reduced. Material that has set inside the mixing tip should not be extruded by force, as this may damage both the cartridge and the mixing tip. Core build-ups which are fabricated using composites, compomers, or resin-modified glasionomer cements must be sealed off before the application of Protemp™ Plus, e.g., using Vaseline.

- First, install a new mixing tip.
- Prior to each application, extrude a small (pea-sized) amount of paste from the mixing tip onto a mixing pad and discard this material.
- Subsequently load the dried impression or vacuum-formed template in the relevant spaces from bottom up.
- Reposition the impression or vacuum-formed template in the mouth.
- The material attains a hard-elastic consistency within 1 min 40 sec after the onset of mixing. The material and the impression or vacuum-formed template must be removed from the mouth within 2 min 50 sec after the onset of mixing.
- Check the setting process by means of the excess material in the mouth.
- Keep the filled mixing tip on the cartridge to serve as a closure until the next use.

Finishing
Option 1 for optimum surface smoothness and gloss
- Completely remove the oxygen inhibition layer using alcohol (e.g., ethanol) once the material is completely cured (no earlier than 5 minutes after the onset of mixing). Polish if desired, e.g., with Sof-Lex™ polishing discs, manu-factured by 3M ESPE. - After the product has cured completely, remove any residues. - If necessary, finish using fine carbide burs and polish, e.g., with Sof-Lex™ polishing discs, manufactured by 3M ESPE. - Do not breathe polishing dust. Use mouth protective device, safety glasses, and aspiration.
Option 2 for faster results
- After removal from mouth and matrix, trim margins using scissors. - Once the material is completely cured (no earlier than 5 minutes after the onset of mixing), completely remove the oxygen inhibition layer using alcohol (e.g., ethanol). Finish using fine carbide burs. - Polish if desired, e.g., with Sof-Lex™ polishing discs, manufactured by 3M ESPE. - Do not breathe polishing dust. Use mouth protective device, safety glasses, and aspiration.

Individualization/Shape Correction
When using a Filtek Flowable Composite and Protemp™ Plus for custom-shaping and shape correction, there usually is no need for a bonding agent. If you intend to use other products, please refer to the pertinent manufacturer's Instructions for Use.
Compatible Products
- Flowable composites, e.g., from the Filtek™ Flowable Composite product line: A1, A2, A3, A3.5, B2, XW (Bleach) - Please comply with the corresponding Instructions for Use. - Protemp™ Plus - Processing: please refer to the information provided under "Application" and "Finishing".

Application of a Filtek Flowable Composite to New Temporary Restorations
In newly manufactured temporary restorations, the Filtek Flowable Composite can be applied either to the finished surface or the unprocessed, but clean, inhibition layer. Bonding is excellent in either case. - Remove any contamination – e.g., from saliva or polishing dust – with alcohol, then dry with a stream of air. Cleaning is not necessary if the inhibition layer is clean or if the temporary restoration is finished and clean. - For reasons of hygiene, dose the Filtek Flowable Composite from the Dispenser onto a mixing pad. - Apply the composite to the temporary restoration in layers of maximally 1 mm using a suitable instrument. - Cure each layer separately with light for 20 sec each. - Trim the composite with carbide burs and polish, e.g., with Sof-Lex™ polishing discs, manufactured by 3M ESPE.

Application of a Filtek Flowable Composite to Temporary Restorations That Were in Use
- Thoroughly roughen the surface of the temporary restoration with rotary instruments. - Remove any contamination – e.g., from saliva or polishing dust – with alcohol, then dry the restoration with a stream of air. - For reasons of hygiene, dose the Filtek Flowable Composite from the Dispenser onto a mixing pad. - Apply the composite to the temporary restoration in layers of maximally 1 mm using a suitable instrument. - Cure each layer separately with light for 20 sec each. - Trim the composite with carbide burs and polish, e.g., with Sof-Lex™ polishing discs, manufactured by 3M ESPE.

Lining of Prefabricated, Temporary Crowns Made of Composite and Metal
Core build-ups which are fabricated using composites, compomers, or resin-modified glasionomer cements must be sealed off before the application of Protemp™ Plus, e.g., using Vaseline. - Shape and cure prefabricated, temporary composite crowns (e.g., Protemp™ Crown, manufactured by 3M ESPE) completely. - Please comply with the corresponding Instructions for Use. - Remove any contamination – e.g., from saliva or polishing dust – with alcohol (e.g., ethanol), then dry with a stream of air. - Apply Protemp™ Plus to the prefabricated, temporary crown; see "Application". - Place the prefabricated crown on the prepared tooth. - The material attains a hard-elastic consistency within 1 min 40 sec after the onset of mixing and must be removed, along with the prefabricated crown, from the mouth within 2 min 50 sec after the onset of mixing. - Check the setting process by observing the residue in the mouth. - Keep the filled mixing tip on the cartridge to serve as a closure until the next use. - Remove the oxygen inhibition layer and finish the residues, see "Finishing".

Cementation
- If the die restoration has been sealed off (e.g., using Vaseline), remove the rest of the sealant thoroughly before the temporary cementation. - Insert the temporary restoration with a commercial temporary cement, e.g., RelyX™ Temp E or RelyX™ Temp NE, both products manufactured by 3M ESPE. - Eugenol-containing cements may impair the setting of composite luting cements to be used at a later time in the final luting step! If a composite cement is to be used for the final luting step, seat the temporary restora-tion with a eugenol-free cement.

After Cementation
Thoroughly examine and explore the sulcus and surrounding dentition. Remove any residuals of the temporary material or cement from the mouth.
Repair
- Provide for mechanical retention by roughening the break surfaces and adjacent areas. - Clean with alcohol and dry with air. - Join together using Protemp™ Plus or a Filtek Flowable Composite. - Cure Protemp™ Plus or the Filtek Flowable Composite as described in the relevant Instructions for Use. - Once the material is fully set, finish as usual.

Cleaning of the Garant Dispenser
- Remove the Garant™ cartridge. - Remove paste that has not set using an alcohol-soaked cloth. - The Dispenser handle and plunger can be autoclaved up to a temperature of 135 °C/275 °F; remove the plunger beforehand. - Glutaraldehyde-based solutions can be used for disinfection.

Storage and Stability
Store the product at 15–25 °C/59–77 °F. Do not store in the refrigerator. Do not use after the expiration date.
Customer Information
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.
Warranty
3M ESPE warrants this product will be free from defects in material and manu-facture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED

WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability
Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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