OptraPol® Next Generation

**Brugsanvisning
Käyttöohjeet
Bruksanvisning
Productinformatie
Οδηγίες Χρήσεως
Kullanım Kılavuzu
Инструкция
Instrukcja stosowania**

**Instructions for Use
Gebrauchsinformation
Mode d'emploi
Istruzioni d'uso
Instrucciones de uso
Instruções de Uso
Bruksanvisning**

- Diamant-poleringsssystem til kompositmaterialer
- Timanttikiillotusjärjestelmä yhdistelmämuovimateriaaleille
- Diamant poleringssystem for kompositmaterialer
- Diamantpolijststelsysteem voor composietmaterialen
- Σύστημα στίλβωσης με διαμάντια για σύνθετα υλικά
- Kompozit materyaller için elmas cilalama sistemi
- Алмазная полировочная система для композитных материалов
- System gumek do polerowania, zawierający rozdrobniony diament

- Diamond polishing system for composite materials
- Diamant-Poliersystem für Compositmaterialien
- Système de polissage pour les matériaux composites
- Sistema di lucidatura diamantato per materiali compositi
- Sistema de pulidores diamantados para materiales composite
- Sistema dotato de partículas de diamante para polimento de compositos
- Diamant-polersystem för kompositmaterial

For dental use only.
Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

Made in Germany
Ivoclar Vivadent AG,
FL-9494 Schaan/Liechtenstein

CE 0123

15.04.2010/WE3

ivoclar
vivadent®
clinical

English

Instructions for Use

Description

OptraPol „Next Generation“ is an effective, high-quality one-step polishing system with diamond particles suitable for polishing all popular direct (e.g. Tetric EvoCeram®) and indirect composite restorative materials to a high lustre. It can also be used for polishing amalgam restorations. The polishers of the OptraPol system are available in the following shapes: small flame, large flame (with red silicone core), cup (with red silicone core) and lens (with red silicone core).

Composition

OptraPol polishing instruments consist of a light grey polishing body and a red silicone core without any polishing effect. The light grey polishing body is composed of silicone rubber and diamond particles. The red silicone core is made of silicone rubber, aluminium oxide, iron oxide and Irgazin red. The polishing body and core are mounted on shanks made of stainless steel.

Indication

OptraPol „Next Generation“ provides effective high-gloss polishing results in just one step. It is indicated for the high gloss polishing of direct and indirect restorations made of all popular types of composite material as well as for amalgam restorations.

Contraindications

If a patient is known to be allergic to any of the above components, the polishers should not be used.

Side effects

None known to date.

Application

The restoration is polished to a high-gloss finish at medium contact pressure (approx. 2N) in a single step. The polishing procedure is performed under water cooling, which also ensures evacuation of polishing residue.

Diamond polishers:

- Rotary speed: 5,000 - 8,000 rpm
- Maximum speed: 10,000 rpm
- Only use with copious water spray (> 50ml/min)
- Use without polishing paste

Proper use

- Use only straight/contra-angle handpieces which are technically and hygienically impeccable and regularly serviced and cleaned.
- The handpiece drives must run smoothly.
- The instruments must be correctly inserted and secured.
- Polishing should be performed with a slight rotary motion to avoid the creation of grooves.
- Jamming/leveraging should be avoided to reduce the risk of instrument fracture.
- Deformed instruments that do not run smoothly should cease to be used.
- From the point of view of safety, protective goggles and breathing protection should always be worn in order to avoid inhaling grinding dust. Efficient dust evacuation must be ensured during the polishing procedure.

Warning

Improper use may lead to poor polishing results and a heightened health risk. This product must only be used by skilled persons.

Rotary speed

The recommended and maximum admissible rotary speeds indicated in the Instructions for Use and on the packaging must be observed. They must not be exceeded. The polishing instruments tend to vibrate if the maximum admissible rotary speed is exceeded. This may lead to the destruction of the polishing instrument, deformation of the shank and/or fracture of the instrument and thus poses a potential danger to the user/practice team and the patient.

Warning

Failure to observe the maximum admissible rotary speed results in a heightened safety risk.

Contact pressure

High contact pressure during polishing will not produce an improved polishing result, but will lead to the increased generation of heat and accelerated wear of the polishing instrument. A mean contact pressure of approx. 2 N is sufficient.

Warning

When polishing intraorally, excessive pressure needs to be avoided in order to prevent overheating and eventually pulp damage. In extreme cases, instrument fracture may result if too much pressure is applied, which can lead to injuries.

Water cooling

In order to avoid unwanted heat generation on the tooth or soft tissue surface, sufficient water cooling (at least 50 ml/min) needs to be ensured. In addition, any residue produced during polishing is rinsed off by the water spray.

Warning

Insufficient water cooling may lead to irreversible damage of the tooth, the restoration and the adjacent tissues.

Storage, disinfection, cleaning and sterilization

- The polishing instruments should be stored in the original packaging at room temperature and protected from dust, moisture, pressure and re-contamination.
- Before using the polishing instruments on the patient for the first time, they must be disinfected and sterilized.
- After having used the polishing instruments in the mouth of the patient, they must be disinfected and cleaned and sterilized immediately.
- When manipulating contaminated instruments, protective clothing and gloves should be worn.
- Mechanically clean the instruments and immerse them in a disinfectant solution.
- Sterilize the polishers in the steam autoclave at 134 °C / 275 °F.
- The products are not suitable for sterilization in the chemi-clave or steam disinfectant.
- Care needs to be taken that the selected disinfectant is suitable for the disinfection of the polishers and that the instruments are disinfected correctly according to the instructions of the manufacturer. Please refer to the instructions of the

manufacturer for details on recommended use, immersion time and suitability of disinfection and cleaning agents. Aggressive substances, excessively long immersion times or incorrect dosage may damage the instruments and reduce the number of re-uses.

Warning

The suitability of the selected method for the disinfection/sterilization of the respective instrument must be verified in advance.

Recommendations to prolong the life of polishing instruments

Mechanical cleaning of the polishing instruments under running tap water and the subsequent sterilization in the autoclave increases the number of times the polishing instruments can be re-used. Alcohol-based disinfectant solutions may reduce reusability.

Storage and shelf life

- Storage temperature 12 - 28 °C/54 - 82 °F
- Expiry date: See imprint on label
-  Sterilizable in the steam sterilizer (autoclave) at the temperature indicated
-  The respective instructions for use can be downloaded from the Ivoclar Vivadent website.

Keep out of the reach of children!
For use in dentistry only!

Date of preparation:

04/2010; Rev. 0

Manufacturer:

Ivoclar Vivadent AG, 9494 Schaan,
Principality of Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.