

RelyX™ Luting Plus Luting Cement



Indications:

For permanent cementation of:

- Porcelain-fused-to-metal crowns and bridges
- Metal crowns, bridges, inlays and onlays (including prefabricated provisional metal crowns).
- All-zirconia or all-alumina strengthened core restorations
- Prefabricated or cast endodontic posts

Cementing orthodontic bands

Tooth Preparation:

- Remove provisional restoration and all temporary cement.
- Clean tooth with oil-free pumice paste, rinse and lightly dry leaving tooth surface moist.
- **DO NOT OVERDRY THE TOOTH**




1

Dispensing and Mixing:

- Remove cap and dispense desired amount of cement onto mixing pad. 1-2 “Clicks” of material are sufficient for a single tooth restoration.
- **Note: To avoid water evaporation and drying out of the cement pastes, dispense and mix cement immediately prior to use.**
- Replace cap securely onto Clicker™ dispenser until an audible click is heard.
- Mix cement with spatula for 20 seconds.


RelyX™ Luting Plus Luting Cement



Load Cement Into Restoration:

- Working time of mixed cement is 2.5 minutes.
- Apply thin layer of cement to the inside surface of the restoration.

2



Placement of Restoration:

- Seat restoration with light pressure.
- Excess cement can be removed after 2 minutes in the mouth with an appropriate instrument (e.g. explorer, scaler).
- Set time is 5 minutes after placement in the mouth.
- Finish restoration and adjust occlusion.

3

3M ESPE

Dental Products, St. Paul, MN 55144-1000 • www.3MESPE.com

3M ESPE Customer Care/MSDS Information:

U.S.A. 1-800-634-2249 and Canada 1-888-363-3685

3M, ESPE, RelyX and Clicker are trademarks of 3M or 3M ESPE AG.

© 3M 2009. All rights reserved.



44-0007-4720-2