



New as of:

11.2006

CEREC 3 / *Scan* / Chairline, inLab

Maintenance Protocol

English



CEREC 3 / *Scan* / Chairline, inLab

Protocol for annual maintenance

| | |
|----------------|----------------|
| Customer | Dealer |
| Address, tel.: | Address, tel.: |

| Serial numbers | Item | Serial no. | | <input checked="" type="checkbox"/> |
|----------------|------|------------|--------------------------|-------------------------------------|
| | 1 | | CEREC 3 acquisition unit | <input type="checkbox"/> |
| | 2 | | CEREC 3 milling unit | <input type="checkbox"/> |
| | 3 | | CEREC <i>Scan</i> | <input type="checkbox"/> |
| | 4 | | inLab | <input type="checkbox"/> |
| | 5 | | CEREC Chairline | <input type="checkbox"/> |
| | 6 | | 3D camera | <input type="checkbox"/> |
| | 7 | | SIROCAM 3 camera | <input type="checkbox"/> |
| | 8 | | SIROCAM 3 module | <input type="checkbox"/> |
| | 9 | | SIROCAM 2 camera | <input type="checkbox"/> |
| | 10 | | SIROCAM 2 module | <input type="checkbox"/> |
| | | | SIDEXIS module | <input type="checkbox"/> |

Prior to starting maintenance, a function check was performed and any defects found were disclosed to the practice staff. These defects will NOT be corrected as part of this maintenance job.

Maintenance was performed according to this maintenance protocol using OEM parts from Sirona.

Date of maintenance

Engineer's name

Signature




Acquisition unit

| Component | Work to be performed | Spare part | REF. | Qty |
|---------------------------|---|----------------------|-----------|-----|
| Cover | Check covers and control unit for damage (cracks, parts broken off). | - | - | - |
| Power supply cable | Check for correct seating and damage. | - | - | - |
| Castors (lockable) | Check condition of castors. Check functioning of parking brake. | Castor with lock | 41 82 403 | 2 |
| Castors (nonlockable) | Check condition of castors. | Castor without lock | 41 82 395 | 2 |
| Fan, unit base | Check functioning. The air current must enter the unit from its base. Remove the filter mat first. Check the air current e.g. with a piece of paper. | Auxiliary fan | 58 97 637 | 1 |
| Filter | Open filter cover in unit base and replace filter. | Filter | 58 33 046 | 1 |
| Fan, PC power supply | Check functioning. The air current must flow out of the PC power supply unit. Check the air current e.g. with a piece of paper. | PC power supply | 58 97 629 | 1 |
| Fan, SIROCAM (optional) | Check functioning. The air current must flow to the front out from the SIROCAM drawer. Check the air current e.g. with a piece of paper (SIROCAM 2 only). | Light module, compl. | 58 55 114 | 1 |
| 3D camera | Check camera for physical damage to the plastic housing and the glass prism. | - | - | - |
| | Image quality check (live and recorded images). Acquire a recorded image of calibration set B. Recalibrate camera if necessary. With 3D software: Have the system calculate the 3D model. Check the waviness of its horizontal surfaces and the height scale accuracy (using the Distance window). With COS software: Check for ripple and height scale accuracy in the Section window. | | | |
| 3D camera | Check for ripple and height scale accuracy in the Section window. | - | - | - |
| Heater, 3D camera | Check functioning. The heater must reach a surface temperature of approx. 50°C within max. 5 minutes after the unit is switched on. | Heater | 58 11 190 | 1 |
| Heater, 3D camera | Check distance between glass prism of 3D camera and surface of heater. A distance of 1 - 2 mm should result with the camera set down. If the distance is too small, compensate by shimming with spacer plate. | Spacer plate | 59 24 134 | 1 |
| SIROCAM camera (optional) | Check exit window for damage. Check the camera separating point for secure seating. | - | - | - |

Acquisition unit

| Component | Work to be performed | Spare part | REF. | Qty |
|---------------------------|---|------------------------|-----------|-----|
| SIROCAM camera (optional) | Check the quality of the SIROCAM image in the application used for color fidelity, blooming and streakiness. Note color differences between live and frozen image. | - | - | - |
| Sirocam heater (optional) | Check functioning. The heater must reach a surface temperature of approx. 50°C within max. 5 minutes after the unit is switched on (SIROCAM 2 only). | SIROCAM heater, compl. | 58 32 907 | 1 |
| SIROCAM heater (optional) | Check distance between exit window of SIROCAM camera and surface of heater. A distance of 1 - 2 mm should result with the camera set down. If the distance is too small, compensate by shimming with spacer plate (SIROCAM 2 only). | Spacer plate | 59 24 134 | 1 |
| Trackball | Remove and clean trackball. Clean wiper ring (only NSI trackball, i.e. ball with adjustable ease of action). | - | - | - |
| Buttons / pedals | Check functioning of foot switch pedal and of right and left mouse buttons. | - | - | - |
| LCD monitor | Check display surface of LCD monitor for dirt and grime. If necessary, clean with a special monitor glass cleanser and a soft cloth. Caution: Do not exert strong force on the display surface! | - | - | - |
| LCD monitor | Check Auto Adjust function (see service manual, Chap. 4.6 "Monitor image flickering"). | - | - | - |
| LCD monitor | Check swivel and tilt functions of LCD monitor. | - | - | - |

Acquisition unit

| Component | Work to be performed | Spare part | REF. | Qty |
|-------------------------|--|-----------------------------|------------------------|-----|
| Keyboard | <p>Switch the unit on and start a word processing program (e.g.: Editor).</p> <p>For spacebars with 3 touch points:</p> <p>Press the following keys in succession:</p>  <p>If the character string kezbrd 90 then appears on the screen and the cursor flashes at the start of the next line, all keyboard functions are operating properly.</p> <p>For spacebars with 4 touch points:</p> <p>Press the following keys in succession:</p> <p>USA</p>  <p>GB</p>  <p>If the character string KEZBRD 90-] then appears on the screen and the cursor flashes at the start of the next line, all keyboard functions are operating properly.</p> <p>Quit the word processing program.</p> <p>Check whether there are any cracks in the foil or its edges have started to come undone (peel off).</p> | Keyboard USA Keyboard GB | 58 73 125 59 00 142 | 1 |
| X-ray module (optional) | Check whether the three X-ray LEDs on the keyboard signal a correct condition following switch-on (alternate flashing at first followed by continuous illumination of the green LED). | - | - | - |
| X-ray module (optional) | Attach X-ray sensor and make the unit ready for exposure via the SIDEXIS application to test it. | - | - | - |

Acquisition unit

| Component | Work to be performed | Spare part | REF. | Qty |
|-----------|--|---------------------|-----------|-----|
| PC drawer | Check whether the drawer of the DVD/CD drive strikes the door when the latter is open. Readjust the Velcro strap of the door fastening if necessary. To do this, take out the metal insert on the inside of the door, loosen the self-adhesive Velcro strap and refasten it with the correct length. Replace the metal insert in the door. | - | - | - |
| Software | Check whether the latest version is installed. If necessary, perform an update after consulting with the user. | COS software update | 58 74 115 | 1 |
| | | 3D software update | 59 77 173 | 1 |

CEREC Chairline

| Component | Work to be performed | Spare part | REF. | Qty |
|-------------------|--|------------|-----------|-----|
| 3D camera | Check camera for physical damage to the plastic housing and the glass prism. | - | - | - |
| 3D camera | Image quality check (live and recorded images). Acquire a recorded image of calibration set B. Recalibrate camera if necessary. | - | - | - |
| | With 3D software: Have the system calculate the 3D model. Check the waviness of its horizontal surfaces and the height scale accuracy (using the Distance window). | | | |
| Heater, 3D camera | Check functioning. The heater must reach a surface temperature of approx. 50°C within max. 5 minutes after the unit is switched on. | Heater | 58 11 190 | 1 |
| Foot switch pedal | Check functioning of foot switch pedal. | - | - | - |

Milling unit

| Component | Work to be performed | Spare part | REF. | Qty |
|-----------------------------|--|---------------------------------|-------------------------------------|--------|
| CC controller board | We recommend replacing the lithium battery on the control board as part of maintenance after 3 years, and then every 3 years thereafter. | Lithium battery | 48 91 219 | 1 |
| Housing | Check general condition. Test function of tank door. Remove any residual adhesive foil from the magnets. Reglue the magnetic locks if necessary. | - | - | - |
| Milling chamber door | Check for cracks and leaks . Check functioning of lock. | Protection window | 41 68 428 | 1 |
| LEDs / keys | Check functioning according to operating instructions. | - | - | - |
| Cover plate for electronics | Check whether the button on the board can be actuated via the download key. Replace the cover plate if necessary. | Cover plate with button | 59 19 662 | 1 |
| Software | Accept data from Info milling unit (serial number, milling time of cylindrical diamond, milling time of tapered diamond, total milling time, total scanning time). | - | - | - |
| Air system | Check pump with pressure switch via Service Software. Check the air outlet at the instruments. | - | - | - |
| Fan | Execute the Fan unit software test (see program option "Service -> Media supply"). | - | - | - |
| Water system | Retrofit the pressure cylinder expansion kit if it has not already been installed. Check pump with pressure switch via Service Software. Check pump, cooler, hoses and flanges for leakage. Check whether there is any leakage of fluid at the shafts and plugs of the milling chamber. Water must be emitted at all nozzles of the gearbox and strike all milling instruments at the same height. | Pressure cylinder expansion kit | 58 85 673 | 1 |
| Water system | We recommend cleaning the water paths of the milling unit with DENTACLEAN (refer to Operating Instructions 59 85 747). | DENTACLEAN | 59 88 543 | 1 |
| Tank | Replace complete tank. Also replace O-ring on suction connection. | Tank cover Tank O-ring | 59 42 763 47 04 792 70 03 189 | 1 each |
| Light barriers | Blow out light barriers (remove dirt). Check ring with slit on gearwheel for dirt and clean if necessary. Perform light barrier test with Service Software. | - | - | - |

Milling unit

| Component | Work to be performed | Spare part | REF. | Qty |
|----------------|--|--|-----------|--------|
| Gear wheels | Determine distortion of the gear wheels by checking their backlash during rotation. To do this, remove the housing cover and turn the relevant gear. The corresponding gear wheels must also turn without backlash. Remove front covers and replace gear wheel stops. | Gear wheel stop | 33 09 098 | 4 |
| Block screw | Replace old block screw with new one. | Type VITA | 46 91 916 | 1 |
| Actuators | Test of stepping motors with Service Software. NOTE: Test procedure is relatively loud. There is no danger of damaging the system. | - | - | - |
| Milling motors | Speed and Touch Test with Service Software. | - | - | - |
| Gearbox | Load and mill Cerec 3D, Patient -> Inlay04, Tooth 16 NOTE: As this design may have been deleted by the user, we recommend bringing a disk containing the design to the maintenance site. | - | - | - |
| Gearbox | Calibrate the milling machine. To do this, check the locking position and the locking button to make sure they do not jam or get stuck (replace locking button if necessary). Following the calibration, insert new instruments. Load the test block with Test, Patient -> Pyramid . Mill it out and measure it. The length of the edge must equal 4.97 mm to 5.07 mm. The difference between the edge lengths must be less than 0.05 mm, the height of the pyramid must equal 5.20 to 5.35 mm. | Step Bur 10, Cylinder Pointed Bur, Cone Bur 10 | | 1 each |
| Cleaning | Clean scanner window carefully (if scanner is installed), reference point, milling chamber, milling chamber window. Replace wiper band on milling chamber window. | Wiper band | 18 74 473 | 1 |

Milling unit

| Component | Work to be performed | Spare part | REF. | Qty |
|-----------|---|---|-----------|-----|
| Scanner | <p>Check for mechanical damage (window). Check the scanner mounting and refasten it if necessary.</p> <p>Sensor - test (Scanner - test) with Service Software. Calibrate scanner.</p> <p>With 3D software: Perform a test scan with the "crown framework" model holder using the "WaxUp" design technique. Have the system calculate the model, do not draw any lines, accept the gray restoration proposal. Activate the mesial view (by double-clicking the light green arrow) and then Cut. Measure the maximum cylinder diameter close to the cut edge using the Distance window. The measured values should be 16,0mm +/- 0,1mm both in vertical (occlusocervical) and horizontal (buccolingual) direction.</p> <p>With COS software: Execute ProbeScan with "pot holder". Measure with Section window of software. To do this, call up the Cursor window with "Window -> Cursor". Measure the pot holder from the left to the right edge by entering the position of the first (x) coordinate. The result must lie within the range of 640 -3/+7 pixels.</p> | - | - | - |
| | Complete the maintenance protocol. Hand the unit over to the customer. | Maintenance Protocol CEREC 3 / Scan / Chairline, inLab | 59 19 712 | 1 |

Maintenance set

The maintenance set (REF: 59 19 670) contains the following components:

| Component | REF. | Qty |
|---|-----------|--------|
| Tank, compl. comprising: cover, filter | 59 42 763 | 1 each |
| tank, | 47 04 792 | |
| O-ring | 70 03 189 | |

| Component | REF. | Qty |
|----------------------|-----------|-----|
| Gear wheel stop | 33 09 098 | 4 |
| VITA block screw | 46 91 916 | 1 |
| Step Bur 10 | | 1 |
| Cylinder Pointed Bur | | 1 |
| Cone Bur 10 | | 1 |
| Wiper band | 18 74 473 | 1 |
| Filter | 58 33 046 | 1 |
| Spacer plate | 59 24 134 | 2 |

The following items should be taken along for any repairs/retrofits/cleanings which may prove necessary:

| Component | REF. | Qty |
|---------------------------------|-----------|-----|
| Pressure cylinder expansion kit | 58 85 673 | 1 |
| Locking button | 59 08 921 | 2 |
| Cover plate with button | 59 19 662 | 1 |
| Castor with lock | 41 82 403 | 2 |
| Castor without lock | 41 82 395 | 2 |
| Auxiliary fan | 58 97 637 | 1 |
| Keyboard (USA) | 58 73 125 | 1 |
| Keyboard (GB) | 59 00 142 | |
| Lithium battery | 48 91 219 | 1 |
| Block changing tool | 27 95 995 | 1 |

| Component | REF. | Qty |
|--------------------------------------|-----------|-----|
| Torque wrench (Bur changing tool) | 41 66 265 | 1 |
| Protection window | 41 68 428 | 1 |
| DENTACLEAN | 59 88 543 | 1 |
| COS software update | 58 74 115 | 1 |
| 3D software update | 59 77 173 | 1 |

The following special tools and documents must be taken along on maintenance calls:

| Component | REF. | Qty |
|--|-----------|-----|
| Special tool for locking button | 59 08 947 | 1 |
| Pot holder for ProbeScan | 58 85 293 | 1 |
| Slide gage with 1/100 mm resolution for measuring test block | - | 1 |
| Service manual | 58 35 694 | 1 |
| List of Spare Parts | 58 62 581 | 1 |

We reserve the right to make any alterations which may be required due to technical improvements.

D 3344.105.01.07.02 11.2006

Sprache: englisch
Ä.-Nr.: 108 114

Printed in Germany
Imprimé en Allemagne

Sirona Dental Systems GmbH

Fabrikstraße 31
64625 Bensheim
Germany
www.sirona.com

in the USA:

Sirona Dental Systems LLC
4835 Sirona Drive, Suite 100
Charlotte, NC 28273
USA

in Canada:

Sirona Canada
3250 Ridgeway Drive - Unit 5
Mississauga, Ontario L5L 5Y6
Canada

Order No

59 19 712 D 3344