NPN 02247152

Strawberry


Your practice is *our* inspiration.™

Temporary Gingival Retraction System

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Aluminum Chloride 15% w/w

Non-medicinal ingredients

Kaolin, Patent Blue V, Purified Water, Strawberry flavor.

PROPERTIES

EXPASYL is a paste for temporary gingival retraction that ensures separation of the marginal gingiva and drying of the sulcus.

The main characteristics of EXPASYL are as follows:

- Gentle to the periodontium: specially formulated to avoid damage to the healthy periodontium, preventing gingival recession or bone resorption.
- Optimal ergonomics and appreciable saving of time: gingival retraction is obtained by a single application of EXPASYL into the sulcus. After 1-2 minutes, EXPASYL (easily visible due to its color) is simply removed by an air and water spray and a dry, retracted sulcus is obtained.
- Painless when EXPASYL is used on healthy periodontium.
- EXPASYL has astringent and hemostatic compression and chemical properties: absence of bleeding helps achieve a dry sulcus.
- Average retraction time: 1-2 minutes, according to intraoral conditions.

INDICATIONS

EXPASYL may be used in indications of sulcular retraction such as:

- Dental impressions.
- Seating of temporary and permanent restorations.
- Restoration of Class II and V cavities.

CONTRA-INDICATIONS

- Known allergy to aluminium chloride.
- Presence of periodontal pockets or furcation.

PRECAUTIONS FOR USE

- Use only in one quadrant of the mouth at a time.
- In the case of intra-sulcular preparation, if the preparation reaches the epithelial attachment, it is necessary to perform a surgical or orthodontic crown lengthening procedure to maintain the biological width.
- Verify there is no periodontal and/or furcation involvement.
- Do not introduce the dispensing tip directly into the sulcus.
- Introduce the paste slowly into the sulcus (2 mm/sec).
- As the paste is particularly hydrophilic, avoid contact of EXPASYL with saliva during the procedure.

- Do not swallow.

Note: RINSE THOROUGHLY and verify that EXPASYL is totally removed from the sulcus, as residue may inhibit set of impression materials.

- In case of contact with the eyes, rinse thoroughly with water.
- Keep out of the reach of children.
- Single-use tips. DISPOSE OF the used tips after use.

INSTRUCTIONS FOR USE

Preparation

- Remove the cap from the cartridge and insert the tip.
- Holding applicator in hand, push plunger release mechanism forward with thumb while pulling back on plunger with opposite hand, slowly turning clockwise.
- Insert the tip first into the chamber of the applicator.
- Gently push the plunger forward until it makes contact with the cartridge.
- To extrude the material, gently squeeze the trigger of the applicator, applying pressure evenly to assure continuous flow.

Note: Tip may be bent, if necessary, or curved tip may be used. After tip is in place, bleed cartridge.


figure 1


figure 2


figure 3


figure 4

Application

- Rinse and slightly dry the prepared tooth.
- Slowly inject the paste into the sulcus, without applying pressure. The point of the tip must create a closed space between the tooth and the marginal edge of the gingiva. Do not place the tip on the gingiva. Gradually straighten the axis of the tip (figures 1 & 2).
- When a tooth is not cut (Class II and V conservatory care for example), incline the tip to a nearly horizontal position to create a closed space between the gingiva and the tooth, then put it slowly back in the upright position while injecting the paste.
- Do not tilt the tip away from the long axis of the tooth (figure 3).
- The cord of paste must have a dry and compact appearance; if not, perform a second injection.
- Place a sufficient quantity of product to totally fill the sulcus in order to obtain adequate retraction (figures 4, 5 & 6).

Waiting Time

EXPASYL must be left in place for 1-2 minutes on average; this time depends on two factors:

- The tonicity of the marginal gingiva: 1 minute for a thin one; 2 minutes for a thick and fibrous one.


figure 5

figure 6

figure 7

figure 8

Note: The waiting time is reduced to 30 seconds in the case of epithelial detachment (after rotary curettage or removal with electrosurge for example), the action of the paste is much quicker and more efficient.

Removal

- Remove EXPASYL with an air and water spray with simultaneous aspiration (figure 7).
- Verify that EXPASYL is completely removed to prevent interaction with the impression material (figure 8).
- Do not use any instrument that may damage the tissues.

CLEANING AND STERILIZATION OF THE APPLICATOR

Cleaning

- After each use, remove the capsule and clean the applicator thoroughly under running water to avoid possible corrosion caused by paste.
- Use a disinfecting solution and follow the directions for dilution and use given by the manufacturer.
- Rinse well and dry.

Sterilization

- Sterilization Unit:
134°C maximum, 2.25 bars,
3 minutes minimum
or
121°C, 1.4 bars, 20 minutes

Note: It is important to remove any residual paste on the dispensing applicator prior to sterilization to prevent corrosion due to the aluminum chloride found in Expasyl.

STORAGE

- After each use, remove the capsule from the applicator and re-cap capsule with its original cap.
- Do not store tips in contact with cartridges. Store new tips in their original packaging, away from cartridges (risk of corrosion).
- Do not use past the expiration date indicated on the outer package.
- Store at a temperature < 25 °C in the provided bag for maximum freshness. Do not refrigerate, freeze or store near excess heat.

DISPOSAL

- Do not discard into sewers or streams.
- Do not discard with household waste.
- Dispose of in accordance with all federal, state and local regulations.

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Manufactured by Pierre Rolland for Kerr. Expasyl is a registered trademark and entire property of Pierre Rolland.
Kerr Corporation**
Sybron Canada Limited