

# RELYX TEMP NE

## DESCRIPTION

- RelyX Temp NE has an extremely thin film thickness and therefore is ideal for precision-fit temporary restorations.
- High adhesion to tooth in combination with easy removal of the temporary restoration.
- Eugenol-free, therefore universally applicable, independent on whether subsequent cementations are performed conventionally or adhesively.
- Cementation excess remains in temporary restorations for the most part, therefore laborious cleaning of the tooth or the teeth is no longer necessary.
- Reduction of end hardness possible by adding Vaseline.

## INDICATIONS

- |                                    |              |
|------------------------------------|--------------|
| • Working time (including mixing)  | 2.30 minutes |
| • Setting time (from begin to mix) | 3.30 minutes |

## DIRECTIONS FOR USE

Mix base and catalyst.

## SIMILAR PRODUCTS

Tempbond NE

## PRESENTATION

**Refill:**

1 x Base Paste 36g

1 x Catalyst 16g

1 x Activator

1 x Mixing Pad