

Crown & Bridge Material

SECTION I - PRODUCT / COMPANY IDENTIFICATION

Product Name: ACCESS®CROWNSynonyms: Crown & Bridge materialPart (Item) Number: 360030, 360031, 360032, 360034, 360035 360037, 360038Company Name: Centrix Inc.Address: 770 River RoadCity, State, Zip: Shelton Ct, 06484Emergency Telephone Number:203-929-5582Information Telephone Number:800-235-5862

SECTION II - COMPOSITION / INFORMATION ON INGREDIENTS

Preparation: Mixture of different dimethacrylates, silicate fillers, pigments, catalyst system, stabilizers Description: Cold-curing two-paste-system for temp. crowns, bridges, inlays and onlays, safety cartridges. Hazardous Components:

Name	CAS-No.	Content	Symbol	R-Categories
BIS-GMA	1565-94-2	35-40 %	Xi	36-38
UDMA	41137-60-4	15-20 %	Xi	36-38

SECTION III - POSSIBLE RISKS

Designation of Risk: Xi, irritant Special Details of Risk (R-categories): 36-38 Detrimental effects/symptoms: Irritant in case of eye, skin and inhalation contact

SECTION IV - FIRST AID MEASURES

Inhalation: Fresh air, in case of indisposition or nausea contact a doctor Skin Contact: Rinse thoroughly with plenty of water Ingestion: Contact immediately a physician Eyes Contact: Rinse thoroughly with water, contact a doctor

SECTION V - FIRE-FIGHTING MEASURES

Suitable Extinguishing Agents: Foam, carbon dioxide, powder fire extinguishers. Unsuitable Extinguishing Agents: Water Unusual Fire and Explosion Hazards: None Special Protective Equipment: None

SECTION VI - ACCIDENTAL RELEASE MEASURES

Personal Precautions: Avoid contact of skin or eyes. Wash hands with plenty of water Environmental Precautions: Do not let the material get in contact with sewerage systems. Dispose in accordance with all federal, state, and local regulations. Methods for Cleaning Up: Take up mechanically, use common absorbents

SECTION VII - HANDLING AND STORAGE

Handling Precautions: Only by qualified staff in dental offices or dental laboratories. Safe Storage Conditions: Store at room temperature in good ventilated rooms and areas Other Precautions: Do not store together with food. Keep out of the reach of children.

SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

Respiratory Protection & Ventilation: Do not inhale vapors or dust Protective Clothing or Equipment: Dental office or laboratory working clothes. Unusual Work/Hygienic Practices: Use of protective gloves and safety glasses is recomended. Avoid prolonged contact skin; wash with plenty of water.

SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Paste / Paste Odor: Dimethacrylate-like odor. Color: Light brown / Yellow Form: Pasty / Pasty Other Data: None.

SECTION X - STABILITY AND REACTIVITY.

Stability: Generally stable. Conditions of Instability: Avoid exposure to elevated temperatures and light. Incompatibility with Various Substances: Radical initiators, e. g. peroxides Hazardous Decomposition Products: HCN, Cyanides Further Particulars: None

SECTION XI - TOXICOLOGICAL INFORMATION.

General Remarks: During long standing product application risks have not become unknown

SECTION XII - ECOLOGICAL INFORMATION.

No ecological damage or impact anticipated upon exposure to environment in quantities as packaged

SECTION XIII - DISPOSAL CONSIDERATIONS.

Waste Disposal: Do not dispose into sink or sewerage systems. Dispose in accordance with all federal, state, and local regulations.

SECTION XIV - TRANSPORT INFORMATION.

Special Provisions for Transport: No hazardous goods withing meaning of transport regulations

SECTION XV - REGULATORY INFORMATION

U.S. FDA Dental Device: Class II product EC Medical Device Directive: Class IIa product

SECTION XVI - OTHER INFORMATION

None.

CAUTION: PRODUCT FOR PROFESSIONAL USE ONLY.

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.

EC



aupistrasse 64 b, 657 19 homein, Germany