




ClearView
SINGLE-USE NASAL HOODS

Nasal Hoods and Scavenging Circuits



Serving the dental industry for more than 45+ years!

**MADE
IN THE
USA**

ClearView™ Classic

Now you can actually see your patient breathing through their nose.

Soft Inner Hood

Creates a desirable nasal/facial seal and is comfortable for patient

Low Profile

Expands clinician's field of view and provides easier access to patient's mouth

Keeps patient from making hood adjustments to see dental team activities, X-rays, intraoral photos, TV monitor, etc.

Hood-in-Hood Design

Minimizes ambient nitrous oxide

Universal Fit

Attaches to most scavenging circuits



Clear Outer Hood

Allows clinician to see fit /seal of soft inner hood

Provides indication that patient is breathing through their nose - slight fogging appears within the hood

Single-Use

Minimizes cross-contamination; saves staff time reprocessing

Latex-Free

Protects patients and staff who have latex allergies

ClearView™ Classic

Order No.	Product Description
Pedo	
33037-9	Unscented Grey
33037-10	Outlaw Orange
33037-11	Sassy Strawberry
33037-12	Birthday Bubblegum
33037-17	Groovy Grape
Adult	
33035-9	Unscented Grey
33035-10	Outlaw Orange
33035-12	Birthday Bubblegum
33035-14	French Vanilla
33035-16	Fresh Mint
Large Adult	
33034-9	Unscented Grey
33034-10	Outlaw Orange
33034-14	French Vanilla
33034-16	Fresh Mint

Broad Selection of Colors, Scents & Sizes

Delights patients of all ages and engages them in the selection process; unscented for those who are scent sensitive.



**Free
Nasal Hood
Samples!**

www.accutron-inc.com
800.531.2221

ClearView™ CO₂ Capnography

All the same features of ClearView™ Classic Nasal Hoods and more!

Angled Side Tube Connection

Tube with Luer Lock port leading to capnograph monitor sample line is angled off to the side, creating less obstruction for clinicians and minimizing the potential for kinks

Luer Lock Port

Allows for easy and secure connection to capnograph monitor sample line



Regulatory Requirements

AAOMS Office Anesthesia Evaluations require the use of a capnograph machine to monitor end-tidal CO₂ when administering moderate sedation, deep sedation and general anesthesia, unless precluded or invalidated by the nature of the patient, procedure or equipment.



ClearView™ CO₂ Capnography

Order No. Product Description

Pedo

43037-9	Unscented Grey
43037-12	Birthday Bubblegum

Adult

43035-9	Unscented Grey
43035-14	French Vanilla

Large Adult


43034-9	Unscented Grey
43034-16	Fresh Mint

ClearView™ Nasal Hoods connect effortlessly to your current scavenging circuit, the guide below will help you determine if an adaptor is necessary...

ClearView™ Adaptor Easy Order Guide

Accutron's ClearView™ Nasal Hoods connect effortlessly to:

- ClearView™ Scavenging Circuits I (43005 & 43006)
- ClearView™ Scavenging Circuit II (43003 & 43004)
- Any existing PIP+™ Scavenging Circuit
- Any Porter Brown Scavenging Circuit (no adaptor required)
- Matrix Scavenging System (has one smooth tube on each side & cone with two connections) requires adaptor part number 35360



If you need assistance identifying your scavenging circuit brand call us, toll-free, at **800.531.2221** We'll be happy to help!

Nitrous Oxide Conscious Sedation



Order No. 43003

ClearView

Scavenging Circuits

- Vacuum control assembly meets NIOSH and ADA proposed flowrate recommendations
- Latex- free components
- X-ray compatible system
- Autoclavable circuit
- White finish complements clinical environment
- Accommodates most flowmeters brands

ClearView™ Scavenging Circuits - for Standard Bag Tee


Order No.	Product Description
43003	ClearView™ Scavenging Circuit II (includes Scavenging Circuit with Vacuum Control Valve, Corrugated Tube, Reservoir Bag and complimentary ClearView Single-Use Nasal Hoods)
43005	ClearView™ Scavenging Circuit I (includes Scavenging Circuit with Vacuum Control Valve and complimentary ClearView Single-Use Nasal Hoods)

All Accutron products are sold through authorized dental distributors.

ClearView™ Scavenging Circuits - for Remote Flow Systems & Accu-Vac

Order No.	Product Description
43004	ClearView™ Scavenging Circuit II (includes Scavenging Circuit Corrugated Tube, Reservoir Bag and complimentary ClearView Single-Use Nasal Hoods)
43006	ClearView™ Scavenging Circuit I (includes Scavenging Circuit and complimentary ClearView Single-Use Nasal Hoods)

ClearView™ Scavenging Circuit Components Diagram (standard bag tee)



Nasal Hoods, Scavenging Circuit & Rubber Goods are Latex Free!

Note: For ClearView™ Scavenging Circuits for RFS and Accu-Vac, parts 25591, 22122 and 25592 are replaced with a 6' Spiral Vacuum Tube (part 27132).