










The scentillating choice for
N₂O delivery... because
patient comfort is
your top priority.



Personal Inhaler Plus+™ Nasal Hoods

-  Reduces risk of cross-contamination
-  Relaxes patients, optimizing operator control
-  Saves valuable staff time
-  8 soothing scents/3 comfortable sizes
-  Made in USA for maximum quality standards
-  Non-latex material minimizes potential for allergic reactions
-  Thermoplastic elastomer reduces concern for nasal hood collapse

Accutron Scavenging Circuit

- ▶ Easy to-read vacuum-flow control gauge
- ▶ Low-profile scavenging assembly expands clinical field of operation
- ▶ Clear exhalation valve hub permits respiratory monitoring
- ▶ All components latex-free
- ▶ Convenient X-ray compatible circuitry
- ▶ Easy-to-clean, autoclavable circuitry
- ▶ Flexible system accommodates all brands of flowmeters
- ▶ Bright white finish provides an updated clinical appearance
- ▶ Made in U.S.A. for maximum quality standards



PIP+™ Nasal Hoods

Accutron's Personal Inhaler Plus+™ single-patient-use nasal hood gives your patients reliable, comfortable, high-quality N₂O delivery.

As the only colored and scented nasal hood manufactured in the USA, Accutron's PIP+ sets the quality standard. The PIP+ is made of thermoplastic elastomer to minimize concern for nasal hood collapse and is certified latex-free, reducing the possibility of allergic reactions.



Packages of 24 brightly colored nasal hoods are offered in eight different scents (plus unscented gray). Available in three sizes—large, medium, small—the PIP+ provides patients with a comfortable fit and a desirable clinical seal to prevent N₂O leakage.

Take special care of your patients, your staff and yourself when you offer your patients the Accutron Personal Inhaler Plus+.

The original flavor-scented nasal hood that pleases patients and eases your procedures.

Order PIP+™ Nasal Hoods in these scents and sizes

Small - 24 Per Package

- 33017-17 Groovy Grape
- 33017-12 Birthday Bubblegum
- 33017-11 Sassy Strawberry
- 33017-10 Outlaw Orange
- 33017-9 Unscented Whisper Gray
- 33017-1 Variety Pack 1 (Orange, Strawberry & Bubblegum)

Nasal hoods are individually wrapped and packaged in 24-count containers. To maintain scent integrity, Variety Pack same-scent hoods are packaged in individual zip-locked bags.

Medium - 24 Per Package

- 33016-16 Fresh Mint
- 33016-15 Perfect Peach
- 33016-14 French Vanilla
- 33016-13 Piña Colada
- 33016-12 Birthday Bubblegum
- 33016-11 Sassy Strawberry
- 33016-10 Outlaw Orange
- 33016-9 Unscented Whisper Gray
- 33016-1 Variety Pack 1 (Orange, Strawberry & Bubblegum)
- 33016-2 Variety Pack 2 (Piña Colada, Vanilla & Orange)
- 33016-3 Variety Pack 3 (Mint, Peach & Vanilla)


Large - 24 Per Package

- 33015-16 Fresh Mint
- 33015-15 Perfect Peach
- 33015-14 French Vanilla
- 33015-9 Unscented Whisper Gray
- 33015-3 Variety Pack 3 (Mint, Peach & Vanilla)

24 nasal hoods per package

See back page for Component Products & Adaptor Easy Order Guide.

Scavenging Circuit I (33010)



Scavenging Circuit

Accutron's **Scavenging Circuit** reduces residual nitrous oxide in treatment rooms and surrounding areas, providing patients and staff with the safest environment possible.

X-ray compatible and autoclavable, the system consists of a clear exhalation valve hub, a vacuum control valve assembly (33010 & 32007 only) and a continuum of tubing. The transparent hub is unique in that it permits visual monitoring of patient respiration during dental procedures. In addition, the hub's low-profile design offers an expanded field of operation, a definite plus since treatment space is often in short supply. The circuit's vacuum control assembly meets NIOSH and ADA proposed scavenging flow rate recommendations. And its gauge is easy


to read and highly visible.

For your
protection.

With its bright white tubing, the circuit creates a fresh clinical appearance as it safely and effectively delivers and removes nitrous oxide. Add an Accutron PIP+ single-patient-use nasal hood to ensure proper facial seal and complete the system.

Accutron Scavenging Circuit—the surest way to minimize patient anxiety and secure a safe and productive work environment for doctor and staff. Latex-free components.

Scavenging Circuit II (32007)



Order No. Description

33010	Scavenging Circuit I: Scavenging Circuit with Vacuum Control Valve and complimentary sampling of single-use PIP+ Nasal Hoods (nasal hoods included are representative of all sizes available)
32007	Scavenging Circuit II: Scavenging Circuit with Vacuum Control Valve, Corrugated Tube, Reservoir Bag and complimentary sampling of single-use PIP+ Nasal Hoods (nasal hoods included are representative of all sizes available)
33008	Scavenging Circuit I for Remote Flow System and Accu-Vac: Scavenging Circuit and complimentary sampling of single-use PIP+ Nasal Hoods (nasal hoods included are representative of all sizes available)
33009	Scavenging Circuit II for Remote Flow System and Accu-Vac: Scavenging Circuit, Corrugated Tube, Reservoir Bag and complimentary sampling of single-use PIP+ Nasal Hoods (nasal hoods included are representative of all sizes available)

See back page for Component Products.

Component Products

Adaptors/Reducers – PIP+™ & Multi-Use Nasal Hoods to Scavengers:

35000	Porter Brown Adaptor
35100	Porter Clean Air Adaptor
35325	Matrx Spiral Tubing Adaptor
35350	Matrx Smooth Tubing Adaptor
35400	McKesson Adaptor
35600	Non-Scavenging System Tube Size Reducing Kit
20305	Side Tube Reducers (2/Pkg)
26250	Cone Ring (Matrx)

Scavenging Circuits with Multi-Use Nasal Hoods:

33011	Scavenger with Large Multi-Use Hood
33012	Scavenger with Medium Multi-Use Hood
33013	Scavenger with Small Multi-Use Hood



Multi-Use Nasal Hoods:


33018	Multi-Use Large Hood
33019	Multi-Use Medium Hood
33020	Multi-Use Small Hood

Scavenging Circuit Components:

21178	Double Connectors for Scavenging System (2/Pkg)
21179	Single Connectors for Nasal Hood (2/Pkg)
21245	Large Y Connector for Scavenging System
21246	Small Y Connector for Scavenging System
22099	Scavenging Hub (clear disk)
22122	Vacuum Gauge
25591	Spiral Vacuum Tube – 3.5'
25592	Spiral Vacuum Tube – 1.5'
25593	5" Large Tube
25594	7 1/4" Small Tube
26253	Tube Slide Double
26254	Single Vacuum Connector
26255	90° Elbow (creates 90° angle off flowmeter for attachment of scavenging circuit)



Scavenging Circuit Parts Diagram



Note: If using RFS & Accu-Vac Scavenging Circuit, parts 25591, 22122 and 25592 are replaced with a 6' Spiral Vacuum Tube (part 26196).

Rubber Goods:

22137	6-foot Corrugated Tube
33014	3-foot Corrugated Tube
33025	Reservoir Bag – 3 liter
33026	Reservoir Bag – 2 liter



Adaptor Easy Order Guide

Porter Brown Scavenging System - 35000

Has 4 tubes - 2 ports in each side, hood within hood bearing the name "Brown."

Porter Clean Air Scavenging System - 35100

Has tube inside tube and wings or vacuum sweeps to draw off excess gas.

Matrx Scavenging System - 35350

Has one smooth tube on each side.

Matrx Scavenging System - 35325

Has spiral, white tubing.

McKesson Scavenging System - 35400

Has clear hood and 1 tube on each side.

Non-Scavenging System Tube Size Reducing Kit - 35600

Reduces tubing to proper size for easy connection.

Accutron Scavenging System does not require an adaptor.

ACCUTRON Inc.
www.accutron-inc.com 1-800-531-2221

Part No. 25558 - Rev. 03

1733 W. Parkside Lane • Phoenix, AZ 85027 • 623-780-2020 • Fax: 623-780-0444